



DGM SI / SA

For Surgical and Anesthetic Instruments

- Single or double door (pass-through) version.
- Economy of the work space: one drying cabinet is enough for loading the volume of up to two washer disinfector capacities.
- Total control for the sterile air into the camera: hepa filtration, overpressure in chamber, indirect UV air treatment during the whole cycle (optional).
- Easy to integrate: lockable doors are reversible and configurable right or left opening, that can be changed during installation.

Model	Capacity	Temperature settings	Time settings, min	Power, kW	Power supply, V/Hz
DGM SI	Up to 18 DIN 1/1 baskets	From ambient to 80 °C	From 1 up to 999 min or continuous	1,7	230/50 other electrical connections on request
DGM SA	Up to 36 anesthesia hoses/18 anesthesia bags and 9 DIN 1/1 baskets			1,75	