



DGM QX 1200

Pre Sterilization Treatment

The DGM-QX-1200 Ultrasonic washer is designed for cleaning surgical tools, laboratory dishware and other medical tools and appliances (in CSSD and local sterilization departments, laboratories, etc.).

- Complete set composition: Ultrasonic cleaner with a sink tap and a spray gun with replaceable nozzles (8 pcs.) and chambers inserts. Hooks for the trays – 2 pcs. Compressor – 1 pcs.
- Perfect decision for express pre sterilization treatment.

Type	Value
External dimensions (D×W×H), mm	600×1400×850
Two chambers, L	60 + 60
Temperature adjustment range, °C	20÷80°C (±3°C), control accuracy 1°C
Operating frequency of the generator, kHz	40 (±5%)
Adjustment range of ultrasounds operation time, min	0÷99 min., control accuracy 1 min.
Power supply, V/Hz	220/50 (other electrical connections on request)

EUROPEAN STANDARD

EN ISO 15883- 1